

Chemilizer Adjustable

Water-Driven, Diaphragm Pump

Description	Product code
20:1 to 64:1	HYDCH9000-AV1
50:1 to 164:1	HYDCH9000-AV2



Maintenance Kits

Chemical Pumps



Change
as needed

Model	Part Number
20:1 to 64:1 (Green)	HYD10095982
50:1 to 164:1 (Blue)	HYD10095989

Check Valve Kit



Change every
6 months

Model	Part Number
20:1 to 64:1	HYD10096924
50:1 to 164:1	HYD10096905

Pump Stem Kit



Change every
6-9 months

Model	Part Number
20:1 to 64:1	HYD10096922
50:1 to 164:1	HYD10096903

Upper Stem O-Ring



Change every
year

Model	Part Number
20:1 to 64:1	HYD10095063
50:1 to 164:1	HYD10095063

Stem O-Ring Gasket



Light Usage:

Change every 6-12 months

Medium Usage:

Change every 6 months

Heavy Usage:

Change every 3-4 months

Model	Part Number
20:1 to 64:1	HYD10098848
50:1 to 164:1	HYD10098633

Pump Rebuild Kit

Change every
6 months



Model	Part Number
20:1 to 64:1	HYD10096920
50:1 to 164:1	HYD10096938

*Images for information purposes only. Kits may have slight variations depending on the model.

**In order for the 50:1 to 1:164 (2.0 to 0.6%) adjustable lower end to work, it must be installed onto a Chemilizer pump that features the brass inserts.

Definitions

Light Usage: Pump is used once a week

Medium Usage: Pump is used 3-4 times each week

Heavy Usage: Pump is used 24 hours a day, seven days a week